

This guide is for all Smartphone users.



How to

yellow.

Use Sensor Modes

The table below details the various AbiBird sensor modes.

ON (Monitoring) To turn the sensor ON for the first time: mode Remove the sensor from its packaging and remove the battery cover at the back of the sensor Push the black button once The sensor will beep to confirm it is on Bluetooth mode Bluetooth mode is a sensor's connecting mode. Use this mode to pair your sensor with your AbiBird app. To activate Bluetooth mode: With the sensor in monitoring mode, press the black button once The sensor will turn blue and flash pink when it detects motion. To deactivate Bluetooth mode: Веер...Веер With the sensor in Bluetooth mode, press the black button once. The blue light will turn off and the sensor will beep twice Alternatively, just wait – the sensor will return to monitoring mode automatically after 5 minutes. OFF (Sleep or Sleep mode is the sensor's OFF mode. Transport) mode Use this mode when you want to turn the sensor OFF, such as when you are transporting it. To activate Sleep mode: With the sensor in monitoring mode, press and hold the black button until the sensor lights up





This guide is for all Smartphone users.



'Call for Help' mode

'Call for Help' mode is the sensor's vibration alarm.

Your cared-for person uses this mode to trigger a 'Call for Help' alert.

To activate 'Call for Help' mode:

 Continuously tap or shake the sensor firmly until it starts beeping and the light turns red.

An alert will be sent to your phone and the sensor's graph will show a 'Call for Help' alert (one or more yellow dots).

